

Quality Control Certificate



Aldosterone RIA (CT)


REF MG13051**LOT** 263710K **2026-11-13**
YYYY-MM-DD

Final Release Results:

		Measured value	
Total Activity	TA	40636	cpm
Maximum Binding/TA	B0/TA	67.94	%

	Concentration	Measured value	Range	B/B0
	pg/mL	cpm	cpm	
CAL 0	0.0	27610	-	100%
CAL 1	29.6	22619	-	82%
CAL 2	58.9	19337	-	70%
CAL 3	107.0	16885	-	61%
CAL 4	312.0	11567	-	42%
CAL 5	611.0	8656	-	31%
CAL 6	1532.0	5745	-	21%

	Measured value	Acceptable range		
	pg/mL	lower limit	target value	upper limit
	pg/mL	pg/mL	pg/mL	pg/mL
CONTROL 1	195.54	135.00	193.00	251.00



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-09-16

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl