

Quality Control Certificate




Erythropoietin (EPO) ELISA

REF NM56011	LOT 3797	2026-06-30 YYYY-MM-DD
--------------------	-----------------	---------------------------------

Final Release Results:

	Concentration mIU/mL	Measured value OD450nm	Measured value OD405nm	OD/ODmax 450nm
CAL A	0	0.000	0.000	0%
CAL B	10.7	0.107	0.031	3%
CAL C	26.1	0.314	0.090	9%
CAL D	51.5	0.653	0.185	19%
CAL E	152	1.860	0.529	54%
CAL F	451	3.474	1.227	100%

	Concentration found mIU/mL	Acceptable range		
		lower limit mIU/mL	target mIU/mL	upper limit mIU/mL
CONTROL 1	19.3	15.3	19.8	24.2
CONTROL 2	177	147	190	232



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-11-26

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl