

Quality Control Certificate



Erythropoietin (EPO) ELISA

REF

NM56011

LOT

3799




2026-06-30

YYYY-MM-DD

Final Release Results:

	Concentration mIU/mL	Measured value OD450nm	Measured value OD405nm	OD/ODmax 450nm
CAL A	0	0.007	0.008	0%
CAL B	10.7	0.163	0.052	5%
CAL C	26.1	0.377	0.115	11%
CAL D	51.5	0.746	0.219	21%
CAL E	152	2.056	0.592	59%
CAL F	451	3.474	1.404	100%

	Concentration found mIU/mL	Acceptable range		
		lower limit mIU/mL	target mIU/mL	upper limit mIU/mL
CONTROL 1	20.0	15.3	19.8	24.2
CONTROL 2	204	147	190	232



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-07-22

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl