

# Quality Control Certificate




## Erythropoietin (EPO) ELISA

<b>REF</b>	<b>NM56011</b>	<b>LOT</b>	<b>3858</b>	<b>2026-12-31</b> YYYY-MM-DD
------------	----------------	------------	-------------	---------------------------------

**Final Release Results:**

	Concentration mIU/mL	Measured value OD450nm	Measured value OD405nm	OD/ODmax 405nm
CAL A	0	0.012	0.009	0%
CAL B	10.7	0.168	0.061	3%
CAL C	26.1	0.412	0.145	8%
CAL D	51.5	0.881	0.301	17%
CAL E	152	2.429	0.812	45%
CAL F	451		1.823	100%

	Concentration found mIU/mL	Acceptable range		
		lower limit mIU/mL	target mIU/mL	upper limit mIU/mL
CONTROL 1	21.6	15.3	19.8	24.2
CONTROL 2	184	147	190	232



**This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.**

Verified by: Dr. Spangenberg

Validation date: 2025-12-30

Accepted

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)