

Quality Control Certificate



Erythropoietin (EPO) ELISA

REF

NM56011

LOT

3890




2027-05-31

YYYY-MM-DD

Final Release Results:

	Concentration mIU/mL	Measured value OD450nm	Measured value OD405nm	OD/ODmax 405nm
CAL A	0	0.012	0.007	0%
CAL B	12.4	0.187	0.068	4%
CAL C	27.0	0.477	0.168	11%
CAL D	50.2	0.906	0.315	20%
CAL E	157	2.543	0.855	54%
CAL F	430		1.590	100%

	Concentration found mIU/mL	Acceptable range		
		lower limit mIU/mL	target mIU/mL	upper limit mIU/mL
CONTROL 1	23.1	17.6	22.0	26.4
CONTROL 2	181	123	174	225



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-05-27

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl