

Quality Control Certificate




Calcitonin ELISA

REF NM59011	LOT 3856	2026-12-31 YYYY-MM-DD
--------------------	-----------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value OD450nm	Measured value OD405nm	OD/ODmax 450nm
CAL A	0	0.020	0.009	1%
CAL B	11.2	0.084	0.029	2%
CAL C	33.4	0.220	0.079	6%
CAL D	112	0.702	0.240	18%
CAL E	329	2.200	0.736	56%
CAL F	1066	3.960	2.491	100%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	23.7	18.1	23.4	28.6
CONTROL 2	419	328	423	518



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-03-12

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl