

Quality Control Certificate



cAMP RIA (Plasma)

REF RE11021

LOT RCA452


2026-11-05
YYYY-MM-DD

Final Release Results:

		Measured value	
Total Activity	TA	61633	cpm
Non-specific binding / TA	NSB/TA	3.2	%
Maximum Binding/TA	(B0-NSB)/TA	33	%

	Concentration pmol/mL	Measured value cpm	Range cpm	(B-NSB) / (B0- NSB)
B0	0	22470	-	100%
CAL A	0.05	21340	-	94%
CAL B	0.15	18814	-	82%
CAL C	0.5	14671	-	62%
CAL D	1.5	9873	-	39%
CAL E	5	5876	-	19%
CAL F	15	3768	-	9%

	Concentration found pmol/mL	Acceptable range		
		lower limit pmol/mL	target pmol/mL	upper limit pmol/mL
CONTROL P	27	19	30	41



Kit Composition:

REF		LOT	
KRCA151	LYO CONTROL P	H387	2027-04
KRCA011	LYO TRACER	H807	2026-11-05
KRCA031	ANTISERUM	H680	2028-03
KRCA111	CAL A-I	G835-Q	2028-08
KRZZ141	CONC SEPREAG	H391	2027-04

REF		LOT	
KRCA281	CONC ASSAYBUF	H692	2028-04
KRCA271	NSB	H681	2028-03
KRCA251	TRIETHYLAMINE	H635	2028-02
KRCA261	ACETANH	H636	2028-02

Valid edition of the instructions for use: V2017-09

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-08-13

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl

Seite 2/2