

Quality Control Certificate



cAMP RIA (Urine)

REF RE11021

LOT RCA453


2026-12-30
YYYY-MM-DD

Final Release Results:

		Measured value	
Total Activity	TA	52339	cpm
Non-specific binding/TA	NSB/TA	1.1	%
Maximum Binding/TA	(B0-NSB)/TA	47	%

	Concentration pmol/mL	Measured value cpm	Range cpm	(B-NSB) / (B0-NSB)
B0	0	25433	-	100%
CAL D	1.5	23218	-	91%
CAL E	5	20118	-	79%
CAL F	15	16051	-	62%
CAL G	50	10665	-	41%
CAL H	150	6589	-	24%
CAL I	450	3525	-	12%

	Concentration found nmol/mL	Acceptable range		
		lower limit nmol/mL	target nmol/mL	upper limit nmol/mL
CONTROL U	2.2	1.8	2.4	3.0



Kit Composition:

REF		LOT	
KRCA171	LYO CONTROL U	H386-J	2027-04
KRCA011	LYO TRACER	H854	2026-12-30
KRCA031	ANTISERUM	H680	2028-03
KRCA111	CAL A-I	H759-A	2028-09
KRZZ141	CONC SEPREAG	H391	2027-04

REF		LOT	
KRCA281	CONC ASSAYBUF	H692	2028-04
KRCA271	NSB	H681	2028-03

Valid edition of the instructions for use: V2017-09

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-10-08

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl