

# Quality Control Certificate



## ARAb Blocking RIA

<b>REF</b> RE21031	<b>LOT</b> RAB280	<b>2026-07-07</b> YYYY-MM-DD
--------------------	-------------------	---------------------------------

### Final Release Results:

			Result		Units
Incubation time	20	3			h
Total Activity (TA)	43868	44945			cpm
Non-specific binding (NSB)	1967	1751			cpm
CONTROL -	8699	8952			cpm
CONTROL +	5984	6078			cpm
NSB/TA	4.5	3.9			% (TA)
CONTROL - / TA	lower limit ≥ 13	upper limit -	20	20	% (TA)
CONTROL +	lower limit 34	upper limit 52	43	43	% Inhibition

### Kit Composition:

<b>REF</b>		<b>LOT</b>	
KRAB013	LYO TRACER	H687A	2026-07-07
KRAB033	NSB REAG	H659-A	2026-10
KRAB153	CONTROL -	H658-A	2026-10
KRAB163	CONTROL +	H660-A	2026-10
KRAB273	CON A SEPH	H567	2026-11
KRAB141	SAMPLEDIL	H519-C	2026-10
KRAR261	BUF	H661-A	2026-07
KRAR253	WASHBUF	H683	2028-03

Valid edition of the instructions for use: V2025-11

<b>This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.</b>	
Verified by:	T. Borchert
Validation date:	2026-04-15
<b>Accepted</b>	
This document has been produced electronically and is valid without a signature.	

