

Quality Control Certificate



ARAb Blocking RIA

REF RE21031

LOT RAB282

2026-10-27
YYYY-MM-DD

Final Release Results:

			Result		Units
Incubation time			20	3	h
Total Activity (TA)			43279	43394	cpm
Non-specific binding (NSB)			1689	1478	cpm
CONTROL -			7946	8650	cpm
CONTROL +			4994	5494	cpm
NSB/TA			3.9	3.4	% (TA)
CONTROL - / TA	lower limit ≥ 13	upper limit -	18	20	% (TA)
CONTROL +	lower limit 34	upper limit 52	46	45	% Inhibition

Kit Composition:

REF		LOT	
KRAB013	LYO TRACER	H803A	2026-10-27
KRAB033	NSB REAG	H659-C	2027-01
KRAB153	CONTROL -	H658-C	2027-01
KRAB163	CONTROL +	H660-C	2027-01
KRAB273	CON A SEPH	H567	2026-11
KRAB141	SAMPLEDIL	H519-E	2027-01
KRAR261	BUF	H743-A	2026-10
KRAR253	WASHBUF	H733	2028-05

Valid edition of the instructions for use: V2025-11

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-08-06

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl