

Quality Control Certificate




Spermatozoa-Ab ELISA

REF RE52029	LOT ESP127	2026-09-30 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration mU/100µL	Measured value OD	manual	
			Range OD	OD/ODmax
CAL A	0	0.151	-	7%
CAL B	31.3	0.331	-	15%
CAL C	62.5	0.571	-	26%
CAL D	125	0.957	-	44%
CAL E	250	1.526	-	70%
CAL F	500	2.184	> 1.300	100%



	Concentration mU/100µL	Acceptable range		
		lower limit mU/100µL	target mU/100µL	upper limit mU/100µL
CONTROL	192	121	186	251

Kit Composition:

REF		LOT	
KESP653	LYO CONTROL	L652	2027-04
KESP515	MTP	L607	2027-03
KESP601	MTP NSB	L608	2027-03
KESP541	CONC ENZCONJ	L630	2027-04
KESP623	LYO CAL	L651	2027-04
KESP753	LYO DILBUF	L615	2027-03
KEZZ701	CONC WASHBUF	L627	2026-09
KEZZ991	TMB SUBS	L541-2	2027-10

REF		LOT	
KEZZ921	TMB STOP	L553	2027-08

Valid edition of the instructions for use: V2015-11 (01)

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2024-11-07

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl