

# Quality Control Certificate




## 17-beta-Estradiol ELISA

<b>REF</b> RE52041	<b>LOT</b> 42K015	<b>2026-12-25</b> YYYY-MM-DD
--------------------	-------------------	---------------------------------

**Final Release Results:**

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL 0	0	2.12	≥1.20	100%
CAL 1	25	1.92	-	91%
CAL 2	100	1.64	-	77%
CAL 3	250	1.30	-	61%
CAL 4	500	1.00	-	47%
CAL 5	1000	0.67	-	32%
CAL 6	2000	0.35	-	17%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL Low	77.1	34.2	62.2	90.2
CONTROL High	262	128	233	338



**This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.**

Verified by: **Dr. Spangenberg**

**Accepted**

Validation date: **2025-02-07**

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)