

# Quality Control Certificate




## Cortisol ELISA

<b>REF</b> RE52061	<b>LOT</b> ECS138	<b>2027-07-31</b> YYYY-MM-DD
--------------------	-------------------	---------------------------------

### Final Release Results:

	Concentration ng/mL	Measured value OD	manual	
			Range OD	OD/ODmax
CAL A	0	2.135	> 1.000	100%
CAL B	20	1.359	-	64%
CAL C	50	0.976	-	46%
CAL D	100	0.640	-	30%
CAL E	200	0.340	-	16%
CAL F	400	0.174	-	8%
CAL G	800	0.090	-	4%



Validation criteria (manual/automated)				
	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	74	44	73	102
CONTROL 2	206	137	196	255

### Kit Composition:

REF		LOT	⌚
KECS511	MTP	K490	2027-07
KECS541	ENZCONJ	K541	2028-02
KECS611	CAL A-G	K538	2027-08
KECS651	CONTROL 1	K539	2027-08
KECS661	CONTROL 2	K540	2027-08
KEZZ991	TMB SUBS	K132-7	2028-09
KIZZ71	TMB STOP	K522	2028-12
KIZZ01	CONC WASHBUF	K459	2029-04

Valid edition of the instructions for use: V2025-08

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-02-02

# Accepted

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)