

# Quality Control Certificate



## LH Urine ELISA

REF RE52111


LOT 06K066-26

2027-11-30  
YYYY-MM-DD

### Final Release Results:

|       | Concentration<br>mIU/mL | Measured<br>value<br>OD | Range<br>OD | OD/ODmax |
|-------|-------------------------|-------------------------|-------------|----------|
| CAL 0 | 0                       | 0.00                    | -           | 0%       |
| CAL 1 | 10                      | 0.16                    | -           | 8%       |
| CAL 2 | 20                      | 0.31                    | -           | 15%      |
| CAL 3 | 40                      | 0.61                    | -           | 30%      |
| CAL 4 | 100                     | 1.23                    | -           | 61%      |
| CAL 5 | 200                     | 2.03                    | ≥ 1.20      | 100%     |

|              | Concentration<br>found<br>mIU/mL | Acceptable range      |                  |                       |
|--------------|----------------------------------|-----------------------|------------------|-----------------------|
|              |                                  | lower limit<br>mIU/mL | target<br>mIU/mL | upper limit<br>mIU/mL |
| Biorad 40462 | 19.9                             | 14.9                  | 27.1             | 39.3                  |
| Biorad 40463 | 42.4                             | 31.2                  | 56.8             | 82.4                  |



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-06-17

**Accepted**

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)