

Quality Control Certificate



FSH ELISA

REF RE52121


LOT 15K125

2027-11-12
YYYY-MM-DD

Final Release Results:

| | Concentration mIU/mL | Measured value OD | Range OD | OD/ODmax |
|-------|-------------------------|-------------------------|-------------|----------|
| CAL 0 | 0 | 0.01 | - | 1% |
| CAL 1 | 5 | 0.15 | - | 9% |
| CAL 2 | 10 | 0.28 | - | 17% |
| CAL 3 | 20 | 0.57 | - | 34% |
| CAL 4 | 50 | 1.20 | - | 72% |
| CAL 5 | 100 | 1.66 | ≥ 1.20 | 100% |

| | Concentration found mIU/mL | Acceptable range | | |
|--------------|----------------------------------|-----------------------|------------------|-----------------------|
| | | lower limit mIU/mL | target mIU/mL | upper limit mIU/mL |
| Biorad 40461 | 5.67 | 2.84 | 5.17 | 7.50 |
| Biorad 40462 | 18.5 | 10.5 | 19.1 | 27.6 |
| Biorad 40463 | 35.0 | 18.2 | 33.2 | 48.1 |



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-01-07

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl