

Quality Control Certificate




Testosterone ELISA

REF RE52151	LOT ETE137	2026-09-30 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration ng/mL	manual		
		Measured value OD	Range	OD/ODmax
CAL A	0	1.784	-	100%
CAL B	0.2	1.574	-	88%
CAL C	0.5	1.375	-	77%
CAL D	1	1.104	-	62%
CAL E	2	0.770	-	43%
CAL F	6	0.339	-	19%
CAL G	16	0.131	-	7%



Validation criteria				
	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target value ng/mL	upper limit ng/mL
CONTROL 1	0.80	0.46	0.77	1.08
CONTROL 2	3.98	2.79	3.99	5.19

Kit Composition:

REF		LOT	
KETE511	MTP	L654	2026-09
KETE541	ENZCONJ	L655	2026-09
KETE611	CAL A	L416-2A	2026-09
KETE611	CAL B-G	L416-2A	2026-09
KETE651	CONTROL 1	L417-2A	2026-09
KETE661	CONTROL 2	L418-2A	2026-09

REF		LOT	
KIZZ01	CONC WASHBUF	L653	2028-04
KEZZ991	TMB SUBS	M840-5	2027-02
KEZZ921	TMB STOP	M920	2026-11

Valid edition of the instructions for use: 2024-01

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: M.Dries

Validation date: 2024-11-20

Accepted

This document has been produced electronically and is valid without a signature.