

Quality Control Certificate




Testosterone ELISA

REF	RE52151	LOT	ETE138	2028-07-31 YYYY-MM-DD
-----	---------	-----	--------	--------------------------

Final Release Results:

	Concentration ng/mL	manual		
		Measured value OD	Range	OD/ODmax
CAL A	0	2.136	-	100%
CAL B	0.2	1.812	-	85%
CAL C	0.5	1.493	-	70%
CAL D	1	1.247	-	58%
CAL E	2	0.868	-	41%
CAL F	6	0.379	-	18%
CAL G	16	0.150	-	7%

Validation criteria				
	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target value ng/mL	upper limit ng/mL
CONTROL 1	0.70	0.49	0.81	1.14
CONTROL 2	4.19	2.91	4.15	5.40



Kit Composition:

REF		LOT	
KETE511	MTP	K938	2028-08
KETE541	ENZCONJ	K939	2028-07
KETE611	CAL A	K941	2028-07
KETE611	CAL B-G	K941	2028-07
KETE651	CONTROL 1	K942	2028-07
KETE661	CONTROL 2	K943	2028-07

REF		LOT	
KIZZZ01	CONC WASHBUF	L862	2028-07
KEZZ991	TMB SUBS	K132-6	2028-09
KEZZ921	TMB STOP	K441	2028-11

Valid edition of the instructions for use: 2024-01

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-09-08

Accepted

This document has been produced electronically and is valid without a signature.