

Quality Control Certificate



DHEA ELISA


REF RE52221**LOT** EDH134 **2028-07-31**
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	manual			e.g.automated/DSX
		Measured value OD	Range* OD	OD/ODmax OD	Measured value OD
CAL A	0.0	1.989	> 1.000	100%	2.316
CAL B	0.37	1.554	-	78%	1.818
CAL C	1.10	1.111	-	56%	1.380
CAL D	3.30	0.696	-	35%	0.894
CAL E	10.0	0.445	-	22%	0.604
CAL F	30.0	0.248	-	12%	0.298

*During automated procedure the optical densities of calibrators and controls may rise due to temperature and specific device conditions. In case of none conformance with the OD ranges stated on the QC Certificate please contact IBL-International or use the validation criteria only!

Validation criteria				
	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	1.8	1.1	1.8	2.5
CONTROL 2	4.6	3.7	4.6	5.6



Kit Composition:

REF		LOT	
KEDH511	MTP	K912	2028-07
KEDH541	ENZCONJ	K945	2028-07
KEDH611	CAL A-F	K925	2028-07
KEDH651	CONTROL 1	K926	2028-07
KEDH661	CONTROL 2	K927	2028-07
KEZZ991	TMB SUBS	K610-6	2029-08
KEZZ921	TMB STOP	K919	2029-09
KIZZ01	CONC WASHBUF	K924	2030-02

REF	LOT	
-----	-----	--

Valid edition of the instructions for use: 2025-03

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-10-08

Accepted

This document has been produced electronically and is valid without a signature.

**IBL International GmbH**
Flughafenstr. 52A, 22335 Hamburg, GermanyTel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl