

Quality Control Certificate




Progesterone ELISA

REF RE52231**LOT** 23K076**2028-01-31**
YYYY-MM-DD

Final Release Results:

	Concentration	Measured value	Range	OD/ODmax
	ng/mL	OD	OD	
CAL 0	0.0	2.65	≥ 1.20	100%
CAL 1	0.30	2.00	-	75%
CAL 2	1.25	1.33	-	50%
CAL 3	2.5	0.98	-	37%
CAL 4	5.0	0.69	-	26%
CAL 5	15.0	0.43	-	16%
CAL 6	40.0	0.31	-	12%

	Concentration found	Acceptable range		
	ng/mL	lower limit	target	upper limit
	ng/mL	ng/mL	ng/mL	ng/mL
Biorad 40461	0.85	0.41	0.75	1.08
Biorad 40462	8.96	4.94	8.97	13.00
Biorad 40463	29.8	16.8		>40



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-09-08

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl