

Quality Control Certificate



Progesterone ELISA

REF RE52231


LOT 23K085

2026-06-30
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL 0	0.0	2.42	≥ 1.20	100%
CAL 1	0.30	1.94	-	80%
CAL 2	1.25	1.25	-	52%
CAL 3	2.5	0.91	-	38%
CAL 4	5.0	0.72	-	30%
CAL 5	15.0	0.38	-	16%
CAL 6	40.0	0.26	-	11%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
Biorad 40441	0.88	0.42	0.76	1.10
Biorad 40442	10.20	4.31	7.86	11.40
Biorad 40443	32.8	21.9		>40



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-08-20

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl