

Quality Control Certificate




Cortisol Urine ELISA

REF	RE52241	LOT	6628A	2027-09-30 YYYY-MM-DD
-----	---------	-----	-------	--------------------------

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL 0	0	2.939	> 1.175	100%
CAL 1	1	2.243	> 0.895	76%
CAL 2	5	1.469	> 0.590	50%
CAL 3	30	0.659	> 0.265	22%
CAL 4	200	0.195	> 0.080	7%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL LOW	3.3	2.7	4.5	6.3
CONTROL HIGH	63.6	44.4	63.4	82.4



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.	
Verified by:	Dr. Spangenberg
Validation date:	2026-07-15
Accepted	
This document has been produced electronically and is valid without a signature.	



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl