

# Quality Control Certificate



## 17-OH-Progesterone Saliva ELISA


**REF****RE52271****LOT****67K066-26****2028-05-22**

YYYY-MM-DD

### Final Release Results:

|       | Calibrator<br>concentration<br>pg/mL | Measured<br>value<br>OD | Range<br>OD | OD/ODmax |
|-------|--------------------------------------|-------------------------|-------------|----------|
| CAL 0 | 0                                    | 1.89                    | $\geq 1.20$ | 100%     |
| CAL 1 | 10                                   | 1.71                    | -           | 90%      |
| CAL 2 | 50                                   | 1.38                    | -           | 73%      |
| CAL 3 | 250                                  | 0.82                    | -           | 43%      |
| CAL 4 | 500                                  | 0.56                    | -           | 30%      |
| CAL 5 | 1000                                 | 0.39                    | -           | 21%      |

|              | Concentration<br>found<br>pg/mL | Acceptable range     |                 |                      |
|--------------|---------------------------------|----------------------|-----------------|----------------------|
|              |                                 | lower limit<br>pg/mL | target<br>pg/mL | upper limit<br>pg/mL |
| CONTROL LOW  | 17.6                            | 11.4                 | 19.0            | 26.6                 |
| CONTROL HIGH | 382                             | 309                  | 515             | 721                  |



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-07-09

## Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.:

E-MAIL:

WEB:

+ 49 (0) 40 532891 -0 Fax: -11

[ibl@tecan.com](mailto:ibl@tecan.com)

[www.tecan.com/ibl](http://www.tecan.com/ibl)