

Quality Control Certificate




Progesterone Saliva ELISA

REF RE52281	LOT EPR142	2026-07-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value OD	manual		e.g. automated EVolyzer
			Range*	OD/ODmax	Measured value OD
CAL A	0	2.082	> 1.000	100%	2.261
CAL B	25	1.639	-	79%	1.872
CAL C	50	1.367	-	66%	1.549
CAL D	250	0.763	-	37%	0.897
CAL E	1000	0.277	-	13%	0.353
CAL F	5000	0.071	-	3%	0.082

*During automated procedure the optical densities of calibrators and controls may rise due to temperature and specific device conditions. In case of none conformance with the OD ranges stated on the QC Certificate please contact IBL-International or use the validation criteria only!



Validation criteria				
	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	57.4	42.5	60.7	78.9
CONTROL 2	576	430	615	799

Kit Composition:

REF		LOT	
KEPR651	CONTROL 1	L551	2026-07
KEPR661	CONTROL 2	L552	2026-07
KEPR511	MTP	L547	2026-07
KEPR541	ENZCONJ	L548	2026-07
KEPR611	CAL A	L549	2026-07
KEPR621	CAL B-F	L550	2026-07
KIZZ01	CONC WASHBUF	L459	2028-01
KEZZ991	TMB SUBS	L541	2027-10

REF		LOT	
KEZZ921	TMB STOP	L458	2027-07

Valid edition of the instructions for use: 2023-10

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2024-09-27

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl