

Quality Control Certificate



Progesterone Saliva ELISA


REF RE52281	LOT EPR143	2026-12-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	manual			e.g. automated EVOlyzer
		Measured value OD	Range* OD	OD/ODmax OD	Measured value OD
CAL A	0	2.177	> 1.000	100%	1.841
CAL B	25	1.694	-	78%	1.488
CAL C	50	1.427	-	66%	1.216
CAL D	250	0.790	-	36%	0.639
CAL E	1000	0.330	-	15%	0.282
CAL F	5000	0.071	-	3%	0.060

*During automated procedure the optical densities of calibrators and controls may rise due to temperature and specific device conditions. In case of none conformance with the OD ranges stated on the QC Certificate please contact IBL-International or use the validation criteria only!

Validation criteria				
	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	62.1	44.8	64.0	83.2
CONTROL 2	708	486	695	903



Kit Composition:

REF		LOT	
KEPR651	CONTROL 1	L827	2026-12
KEPR661	CONTROL 2	L828	2026-12
KEPR511	MTP	L823	2026-12
KEPR541	ENZCONJ	L824	2026-12
KEPR611	CAL A	L825	2026-12
KEPR621	CAL B-F	L826	2026-12
KIZZ01	CONC WASHBUF	L803	2028-06
KEZZ991	TMB SUBS	L541-4	2027-10

REF		LOT	
KEZZ921	TMB STOP	L861	2028-01

Valid edition of the instructions for use: 2024-12

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: M.Dries

Validation date: 2025-02-24

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl