

Quality Control Certificate



Aldosterone ELISA


REF RE52301

LOT 89K066-2

2027-06-30
YYYY-MM-DD

Final Release Results:

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.02	≥ 1.20	100%
CAL B	20	1.54	-	76%
CAL C	80	0.85	-	42%
CAL D	200	0.42	-	21%
CAL E	500	0.17	-	8%
CAL F	1000	0.10	-	5%



	Concentration found pg/mL	Target Value pg/mL	Acceptable range lower limit upper limit pg/mL pg/mL	
Control Low	49.0	50.7	27.9	73.5
Control High	542	542	298	786

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-07-09

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl