

# Quality Control Certificate




## Aldosterone ELISA

<b>REF</b> RE52301	<b>LOT</b> 89K085	<b>2026-07-31</b> YYYY-MM-DD
--------------------	-------------------	---------------------------------

**Final Release Results:**

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	1.76	≥ 1.20	100%
CAL B	20	1.21	-	69%
CAL C	80	0.74	-	42%
CAL D	200	0.43	-	24%
CAL E	500	0.24	-	14%
CAL F	1000	0.15	-	9%



	Concentration found pg/mL	Target Value pg/mL	Acceptable range lower limit upper limit pg/mL pg/mL
Control Low	50.5	52.7	29.0 76.4
Control High	376	442	243 641

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-09-09

**Accepted**

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)