

Quality Control Certificate



Aldosterone Reagent Set Urine

REF	RE52307	LOT	89KU085	2026-08-31 YYYY-MM-DD
-----	---------	-----	---------	--------------------------

Final Release Results:

This Aldosterone Reagent Set Urine is designed for urine pre-treatment before quantification of Aldosterone with RE52307. For instructions for use see IFU RE52301, chapter 5.2.

Dieses Aldosterone Reagent Set Urine dient der Vorbehandlung von Urinproben zur Messung von Aldosteron mit RE52307. Gebrauchsinformationen finden sich in IFU RE52301, Kapitel 5.2.

Component/ Komponente	Release Reagent	Neutralization Buffer	Dilution Buffer
Number x Volume/ Anzahl x Volumen	1 x 3 ml	1 x 3 ml	2 x 25 ml
Lot Number/ Chargenbezeichnung	HC04486357-3	89NB076	89DB076
Expiration Date/ Haltbar bis	31.08.2028	31.07.2028	31.07.2028

BIORAD Lyphochek Quantitative Urine	Mean/Mittelwert (pg/mL)	Acceptance Range/Akzeptanzbereich (pg/mL)	Observed Conc./Gemessene Konz. (pg/mL)
63551*	39387	25602-53172	37734
63552*	108861	70760-146962	108186

* We recommend to dilute the BIORAD Lyphochek Control after preparation 1:10/
Wir empfehlen die BIORAD Lyphochek Kontrolle nach der Vorbereitung 1:10 zu verdünnen.

This product has been tested successfully by the Quality Control Department and was released for sale according

Verified by: Dr. Spangenberg

Validation date: 2025-08-13

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl