

Quality Control Certificate



Cortisol Saliva ELISA


REF RE52611**LOT** ECO172**2028-02-29**

YYYY-MM-DD

Final Release Results:

	Concentration $\mu\text{g/dL}$	Measured value OD	manual	
			Range OD	OD/ODmax
CAL A	0	2.437	> 1.000	100%
CAL B	0.015	2.136	-	88%
CAL C	0.04	1.738	-	71%
CAL D	0.17	1.146	-	47%
CAL E	0.70	0.689	-	28%
CAL F	3.00	0.289	-	12%

	Concentration found $\mu\text{g/dL}$	Acceptable range		
		lower limit $\mu\text{g/dL}$	target $\mu\text{g/dL}$	upper limit $\mu\text{g/dL}$
CONTROL 1	0.066	0.040	0.081	0.121
CONTROL 2	0.55	0.43	0.61	0.79



Kit Composition:

REF		LOT	
KECO651	CONTROL 1	K754-2A	2028-03
KECO661	CONTROL 2	K755-2A	2028-03
KECO511	MTP	K906	2028-02
KECO541	ENZCONJ	K920	2028-09
KECO611	CAL A	K752A-2A	2028-03
KECO621	CAL B-F	K753-2A	2028-03
KIZZ01	CONC WASHBUF	K944	2030-02

REF		LOT	
KEZZ991	TMB SUBS	K610-4	2029-08
KEZZ921	TMB STOP	K903	2029-08

Valid edition of the instructions for use: 2024-02

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-09-23

Accepted

This document has been produced electronically and is valid without a signature.

**IBL International GmbH**
Flughafenstr. 52A, 22335 Hamburg, GermanyTel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl