

Quality Control Certificate




Testosterone Saliva ELISA

REF RE52631	LOT EST155	2027-07-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.104	-	100%
CAL B	10	1.717	-	82%
CAL C	30	1.497	-	71%
CAL D	100	1.046	-	50%
CAL E	300	0.555	-	26%
CAL F	900	0.226	-	11%



	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	25.8	11.8	29.5	47.3
CONTROL 2	209	176	221	265

Kit Composition:

REF	LOT		REF	LOT	
30119833	MTP	K259 2027-07	KIZZZ01	CONC WASHBUF	K272 2029-01
KEST541	CONC ENZCONJ	K276 2027-07	KEZZ991	TMB SUBS	K132-2 2028-09
KEST611	CAL A	L565A-2A 2027-07	KEZZ921	TMB STOP	K270 2028-07
KEST621	CAL B-F	L565A-2A 2027-07			
KEST651	CONTROL 1	L566A-2A 2027-07			
KEST661	CONTROL 2	L567A-2A 2027-07			

Valid edition of the instructions for use: 2025-06

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2025-09-17

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl