

Quality Control Certificate




Testosterone Saliva ELISA

REF RE52631**LOT** EST156 **2028-05-31**
YYYY-MM-DD

Final Release Results:

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.246	-	100%
CAL B	10	1.903	-	85%
CAL C	30	1.585	-	71%
CAL D	100	1.067	-	48%
CAL E	300	0.588	-	26%
CAL F	900	0.264	-	12%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	29.3	11.4	28.4	45.5
CONTROL 2	203	167	209	250



Kit Composition:

REF		LOT	
30119833	MTP	K793	2028-05
KEST541	CONC ENZCONJ	K791	2028-05
KEST611	CAL A	K788	2028-05
KEST621	CAL B-F	K788	2028-05
KEST651	CONTROL 1	K789	2028-05
KEST661	CONTROL 2	K790	2028-05

REF		LOT	
KIZZZ01	CONC WASHBUF	K787	2029-11
KEZZ991	TMB SUBS	K610-3	2029-08
KEZZ921	TMB STOP	K515	2028-12

Valid edition of the instructions for use: 2025-06

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-07-03

Accepted

This document has been produced electronically and is valid without a signature.

IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, GermanyTel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl