

# Quality Control Certificate



## DHEA Saliva ELISA


<b>REF</b> RE52651	<b>LOT</b> EDE131	<b>2027-03-31</b> YYYY-MM-DD
--------------------	-------------------	---------------------------------

**Final Release Results:**

	Concentration pg/mL	manual		
		Measured value OD	Range* OD	OD/ODmax
CAL A	0	1.925	> 1.000	100%
CAL B	20	1.632	-	85%
CAL C	60	1.223	-	64%
CAL D	180	0.929	-	48%
CAL E	540	0.594	-	31%
CAL F	2160	0.278	-	14%

\*During automated procedure the optical densities of calibrators and controls may rise due to temperature and specific device conditions. In case of none conformance with the OD ranges stated on the QC Certificate please contact IBL-International or use the validation criteria only!

Validation criteria (manual/automated)				
	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	94	61	102	143
CONTROL 2	383	292	417	542



**Kit Composition:**

REF		LOT	
KEDE511	MTP	L962	2027-03
KEDE541	ENZCONJ	L984	2027-04
KEDE611	CAL A	L975	2027-04
KEDE621	CAL B-F	L975	2027-04
KEDE651	CONTROL 1	L976	2027-04
KEDE661	CONTROL 2	L977	2027-04
KEZZ991	TMB SUBS	L919	2028-06
KEZZ921	TMB STOP	L861	2028-01
KIZZ01	CONC WASHBUF	L998	2028-10

REF	LOT	

Valid edition of the instructions for use: 2024-08

**This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.**

Verified by: Dr. Spangenberg

**Accepted**

Validation date: 2025-05-16

This document has been produced electronically and is valid without a signature.