

Quality Control Certificate



Estrone Saliva ELISA

REF RE52681


LOT 98K076

2027-07-31
YYYY-MM-DD

Final Release Results:

	Calibrator concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.56	≥ 1.20	100%
CAL B	3	2.42	-	95%
CAL C	12.3	2.16	-	84%
CAL D	37	1.65	-	64%
CAL E	111	1.10	-	43%
CAL F	333	0.74	-	29%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
Control 1	7.5	5.08	9.2	13.4
Control 2	171	113	174	235



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-07-23

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.:
E-MAIL:
WEB:

+ 49 (0) 40 532891 -0 Fax: -11
ibl@tecan.com
www.tecan.com/ibl