

Quality Control Certificate



Chromogranin A ELISA

REF RE53071


LOT AA541

2027-09-30
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.053	-	2%
CAL B	25.6	0.123	-	5%
CAL C	101.0	0.290	-	11%
CAL D	260.5	0.614	-	24%
CAL E	795.8	1.489	-	58%
CAL F	2166.8	2.560	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	31.9	20.0	28.6	37.1
CONTROL 2	154.1	104.0	148.6	193.1




This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-12-10

Accepted

This document has been produced electronically and is valid without a signature.

 IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl