

Quality Control Certificate



Chromogranin A ELISA

REF RE53071


LOT AA946

2028-02-29
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.063	-	2%
CAL B	29.0	0.232	-	8%
CAL C	95.6	0.672	-	23%
CAL D	241.7	1.104	-	37%
CAL E	964.8	2.542	-	85%
CAL F	2025.5	2.983	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	38.7	27.5	39.3	51.1
CONTROL 2	107.9	78.6	112.3	146.0




This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-03-17

Accepted

This document has been produced electronically and is valid without a signature.

 IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl