

Quality Control Certificate



Chromogranin A ELISA

REF RE53071


LOT Y309

2027-03-31
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.088	-	4%
CAL B	26.5	0.200	-	9%
CAL C	133.8	0.423	-	20%
CAL D	221.9	0.534	-	25%
CAL E	782.1	1.115	-	53%
CAL F	2117.7	2.114	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	37.0	28.2	40.3	52.4
CONTROL 2	154.8	101.8	145.4	189.0




This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-04-10

Accepted

This document has been produced electronically and is valid without a signature.

 IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl