

Quality Control Certificate



Osteocalcin (1-43/49) specific ELISA

REF RE53141


LOT Y269

2027-03-31
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL 1	0	0.061	-	2%
CAL 2	5.6	0.236	-	9%
CAL 3	13.6	0.468	-	19%
CAL 4	20.3	0.698	-	28%
CAL 5	40.0	1.535	-	62%
CAL 6	76.3	2.495	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	8.5	5.8	7.7	9.6
CONTROL 2	33.7	22.5	30.0	37.5




This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-06-18

Accepted

This document has been produced electronically and is valid without a signature.

 IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl