

# Quality Control Certificate



## Insulin ELISA

REF RE53171


LOT 44K015

2026-05-31  
YYYY-MM-DD

### Final Release Results:

	Concentration μIU/mL	Measured value OD	Range OD	OD/ODmax
CAL 0	0	0.06	-	3%
CAL 1	6.25	0.12	-	5%
CAL 2	12.5	0.22	-	9%
CAL 3	25	0.44	-	19%
CAL 4	50	1.00	-	43%
CAL 5	100	2.33	≥ 1.20	100%

	Concentration found μIU/mL	Acceptable range		
		lower limit μIU/mL	target μIU/mL	upper limit μIU/mL
Biorad 40441	17.0	9.10	16.6	24.0
Biorad 40442	71.7	41.2		> 100



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-04-01

# Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)