

Quality Control Certificate




Melatonin ELISA

REF RE54021	LOT EME181	2027-07-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.574	> 1.000	100%
CAL B	3.0	2.292	-	89%
CAL C	9.0	2.028	-	79%
CAL D	27	1.316	-	51%
CAL E	82	0.584	-	23%
CAL F	240	0.211	< 0.500	8%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	7.9	4.8	8.8	12.7
CONTROL 2	57.1	43.4	58.7	74.0



Kit Composition:

REF		LOT	
KEME651	LYO CONTROL 1	K196	2027-07
KEME661	LYO CONTROL 2	K197	2027-07
KEME541	CONC ENZCONJ	K194	2027-08
KEME611	LYO CAL A-F	K195	2027-07
KEZZ811	MTP	K191	2028-06
KEME535	LYO BIOTIN	K192	2027-07

REF		LOT	
KEME523	LYO ANTISERUM	K193	2027-07
KIZZ01	CONC WASHBUF	K198	2028-12
KEZZ961	PNPP SUBS	K180	2027-09
KEZZ951	PNPP STOP	K189	2028-05
KEME761	EXTRCOL	L640	

Valid edition of the instructions for use: V2024-11

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2025-09-26

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl