

Quality Control Certificate



Melatonin ELISA

REF RE54021


LOT EME182

2028-04-30
YYYY-MM-DD

Final Release Results:

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.805	> 1.000	100%
CAL B	3.0	2.496	-	89%
CAL C	8.0	2.214	-	79%
CAL D	22	1.498	-	53%
CAL E	64	0.671	-	24%
CAL F	245	0.258	< 0.500	9%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	7.5	4.7	8.5	12.4
CONTROL 2	44.1	33.6	45.4	57.2



Kit Composition:

REF		LOT	
KEME651	LYO CONTROL 1	K728	2028-05
KEME661	LYO CONTROL 2	K729	2028-05
KEME541	CONC ENZCONJ	K726	2028-05
KEME611	LYO CAL A-F	K727	2028-05
KEZZ811	MTP	K191	2028-06
KEME535	LYO BIOTIN	K724	2028-04

REF		LOT	
KEME523	LYO ANTISERUM	K725	2028-04
KIZZ01	CONC WASHBUF	K731	2029-10
KEZZ961	PNPP SUBS	K732	2028-05
KEZZ951	PNPP STOP	K189	2028-05
KEME761	EXTRCOL	K730	

Valid edition of the instructions for use: 2024-11

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Accepted

Validation date: 2026-06-29

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl