

Quality Control Certificate




Melatonin-Sulfate Urine ELISA

REF RE54031	LOT EMS207	2027-01-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0.0	2.551	-	100%
CAL B	2.0	2.041	-	80%
CAL C	6.0	1.341	-	53%
CAL D	18.0	0.657	-	26%
CAL E	54.0	0.246	-	10%
CAL F	162	0.078	-	3%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	5.4	3.7	5.4	7.0
CONTROL 2	25	18	25	33



Kit Composition:

REF		LOT	
KEMS651	CONTROL 1	K315	2027-01
KEMS661	CONTROL 2	K316	2027-01
KEZZ811	MTP	K279	2028-07
KEMS521	ANTISERUM	K313	2027-02
KEMS541	CONC ENZCONJ	K321	2027-02
KEMS611	CAL A-F	K314	2027-01
KEMS751	ASSAYBUF	K306	2027-07

REF		LOT	
KEZZ701	CONC WASHBUF	K280	2027-07
KEZZ991	TMB SUBS	K132-2	2028-09
KEZZ921	TMB STOP	L132	2027-02

Valid edition of the instructions for use: 2023-03

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2025-10-01

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl