

Quality Control Certificate



Melatonin-Sulfate Urine ELISA

REF RE54031


LOT EMS208

2027-11-30
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0.0	2.301	-	100%
CAL B	2.0	1.948	-	85%
CAL C	6.0	1.468	-	64%
CAL D	18.0	0.714	-	31%
CAL E	54.0	0.302	-	13%
CAL F	162	0.093	-	4%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	5.1	3.7	5.3	7.0
CONTROL 2	23	16	22	29



Kit Composition:

REF		LOT	
KEMS651	CONTROL 1	K801	2027-11
KEMS661	CONTROL 2	K802	2027-11
KEZZ811	MTP	K279	2028-07
KEMS521	ANTISERUM	K798	2027-12
KEMS541	CONC ENZCONJ	K799	2028-01
KEMS611	CAL A-F	K800	2027-11
KEMS751	ASSAYBUF	K803	2028-05

REF		LOT	
KEZZ701	CONC WASHBUF	K804	2028-05
KEZZ991	TMB SUBS	K132-7	2028-09
KEZZ921	TMB STOP	K441	2028-11

Valid edition of the instructions for use: 2023-03

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-07-20

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl