

# Quality Control Certificate



## Melatonin direct Saliva ELISA

REF RE54041


LOT ESM147

2027-02-28  
YYYY-MM-DD

### Final Release Results:

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.076	> 1.000	100%
CAL B	0.5	1.883		91%
CAL C	1.5	1.538		74%
CAL D	5.0	1.219		59%
CAL E	15.0	0.650		31%
CAL F	50.0	0.438		21%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	2.9	1.3	2.6	3.8
CONTROL 2	35.6	20.3	31.2	42.2



### Kit Composition:

REF		LOT		REF		LOT	
KEZZ811	MTP	K305	2028-07	KESM535	BIOTIN	K274	2027-02
KESM611	CAL A	K308	2027-02	KESM541	ENZCONJ	K318	2027-02
KESM621	CAL B-F	K309	2027-02	KEZZ751	CONC WASHBUF	K307	2027-08
KESM651	CONTROL 1	K310	2027-02	KEZZ991	TMB SUBS	K132-4	2028-09
KESM661	CONTROL 2	K311	2027-02	KEZZ921	TMB STOP	K270	2028-07
KESM521	ANTISERUM	K275	2027-02				

Valid edition of the instructions for use: 2023-11

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-10-23

# Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)