

Quality Control Certificate



FT3 ELISA

REF **RE55231**


LOT **6389A**

2026-09-30
YYYY-MM-DD

Final Release Results:

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0.0	3.035	> 1.215	100%
CAL B	1.6	2.560	> 1.025	84%
CAL C	3.3	1.971	> 0.790	65%
CAL D	6.0	1.123	> 0.450	37%
CAL E	9.5	0.678	> 0.270	22%
CAL F	19.0	0.280	> 0.110	9%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
BIORAD 40461	2.51	1.74	2.49	3.23
BIORAD 40462	4.89	3.37	4.81	6.25
BIORAD 40463	8.10	4.98	7.12	9.25



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-03-28

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11

E-MAIL: ibl@tecan.com

WEB: www.tecan.com/ibl