

Quality Control Certificate



FT4 ELISA

REF RE55241


LOT 6611A

2028-01-31
YYYY-MM-DD

Final Release Results:

	Concentration ng/dL	Measured value OD	Range OD	OD/ODmax
CAL A	0.00	2.132	> 0.855	100%
CAL B	0.47	1.782	> 0.715	84%
CAL C	1.30	1.320	> 0.530	62%
CAL D	2.40	0.720	> 0.290	34%
CAL E	4.60	0.307	> 0.125	14%
CAL F	7.60	0.111	> 0.045	5%

	Concentration found ng/dL	Acceptable range		
		lower limit ng/dL	target ng/dL	upper limit ng/dL
BIORAD 40471	0.90	0.67	0.96	1.25
BIORAD 40472	1.94	1.34	1.92	2.50
BIORAD 40473	3.28	2.20	3.15	4.09



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-08-27

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl