

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Chlamydia pneumoniae IgG**
 Prod.No.: **RE57041**
 Lot: **CHLPG-201**
 Expiry: **2027-02-28**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0510-201	2025-08	2027-02-28
CONJ	1	KG0510-201	2025-08	2027-02-28
CONTROL -	1	NG0510-201	2025-08	2027-02-28
CUT OFF	1	CG0510-201	2025-08	2027-02-28
CONTROL +	1	PG0510-201	2025-08	2027-02-28
WASH BUF 20x	1	W-1118	2025-08	2027-08-31
SOLN STOP	1	S-1049	2025-07	2027-07-31
SUB TMB	1	T-078	2025-04	2028-04-28
DIL	1	D-493	2025-05	2027-05-31

Component	OD Value	OD Validation Criteria
CONTROL -	0.019	< 0.200 < Cut-off
CUT OFF	0.463	0.150-1.300
CONTROL +	1.014	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	19.92	10.96 - 28.88

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-09-10

Approved by: 
 Quality Control

Approved on: 2025-09-10



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany