

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Chlamydia pneumoniae IgM**

Prod.No.: **RE57051**

Lot: **CHLPM-202**

Expiry: **2027-04-30**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MM0510-202	2025-10	2027-04-30
<b>CONJ</b>	1	KM0510-202	2025-10	2027-04-30
<b>CONTROL -</b>	1	NM0510-202	2025-10	2027-04-30
<b>CUT OFF</b>	1	CM0510-202	2025-10	2027-04-30
<b>CONTROL +</b>	1	PM0510-202	2025-10	2027-04-30
<b>WASH BUF 20x</b>	1	W-1125	2025-10	2027-10-31
<b>SOLN STOP</b>	1	S-1056	2025-09	2027-09-30
<b>SUB TMB</b>	1	T-078	2025-04	2028-04-28
<b>DIL M</b>	1	DM-348	2025-10	2027-10-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.017	< 0.200 < Cut-off
<b>CUT OFF</b>	0.583	0.150-1.300
<b>CONTROL +</b>	1.188	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	20.71	11.39 - 30.02

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-11-21

Approved by: \_\_\_\_\_

Quality Control

Approved on: 2025-11-21



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany