

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Rubella Virus IgG**  
Prod.No.: **RE57081**  
Lot: **RUBG-081**  
Expiry: **2027-08-31**  
IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0400-081	2026-02	2027-08-31
<b>CONJ</b>	1	KG0400-081	2026-02	2027-08-31
<b>CAL A</b>	1	SAG0400-081	2026-02	2027-08-31
<b>CAL B</b>	1	SBG0400-081	2026-02	2027-08-31
<b>CAL C</b>	1	SCG0400-081	2026-02	2027-08-31
<b>CAL D</b>	1	SDG0400-081	2026-02	2027-08-31
<b>WASH BUF 20x</b>	1	W-1137	2026-02	2028-02-29
<b>SOLN STOP</b>	1	S-1069	2026-02	2028-02-29
<b>SUB TMB</b>	1	T-079	2025-12	2028-12-10
<b>DIL</b>	1	D-499	2026-01	2028-01-31

Component	Concentration	OD Value	OD Validation Criteria
<b>CAL A</b>	0 IU/mL	0.009	< 0.200
<b>CAL B</b>	10 IU/mL	0.628	> 0.200 & > Standard A
<b>CAL C</b>	50 IU/mL	1.306	> 0.700 & > Standard B
<b>CAL D</b>	100 IU/mL	1.699	> 1.100 & > Standard C

Final Report issued: 2026-03-18

Approved by: \_\_\_\_\_

Quality Control

Approved on: 2026-03-18



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany