



**GOLD
STANDARD
DIAGNOSTICS**



Certificate of Analysis (CoA)

Product: **Toxoplasma gondii IgG**

Prod.No.: **RE57101**

Lot: **TOXG-090N**

Expiry: **2026-07-31**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0460-090	2025-01	2026-07-31
CONJ	1	KG0460-090	2025-01	2026-07-31
CAL A	1	SAG0460-090	2025-01	2026-07-31
CAL B	1	SBG0460-090	2025-01	2026-07-31
CAL C	1	SCG0460-090	2025-01	2026-07-31
CAL D	1	SDG0460-090	2025-01	2026-07-31
WASH BUF 20x	1	W-1102	2024-10	2026-10-31
SOLN STOP	1	S-1033	2024-09	2026-09-30
SUB TMB	1	T-077	2022-09	2027-03-31
DIL	1	D-488	2024-08	2026-08-31

Component	Concentration	OD Value	OD Validation Criteria
CAL A	0 IU/mL	0.019	< 0.200
CAL B	50 IU/mL	0.939	> 0.300 & > Standard A
CAL C	100 IU/mL	1.608	> 0.500 & > Standard B
CAL D	200 IU/mL	2.431	> 1.000 & > Standard C

Final Report issued: 2025-01-24

Approved by: *Hoe*
Quality Control



IBL International GmbH
Flughafenstrasse 52a, 22335 Hamburg
Germany

Approved on: 2025-01-24



Gold Standard Diagnostics Frankfurt GmbH
Waldstrasse 23 A6
63128 Dietzenbach, Germany