

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Toxoplasma gondii IgG**

Prod.No.: **RE57101**

Lot: **TOXG-092**

Expiry: **2026-10-31**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0460-092	2025-04	2026-10-31
<b>CONJ</b>	1	KG0460-092	2025-04	2026-10-31
<b>CAL A</b>	1	SAG0460-092	2025-04	2026-10-31
<b>CAL B</b>	1	SBG0460-092	2025-04	2026-10-31
<b>CAL C</b>	1	SCG0460-092	2025-04	2026-10-31
<b>CAL D</b>	1	SDG0460-092	2025-04	2026-10-31
<b>WASH BUF 20x</b>	1	W-1108	2025-02	2027-02-28
<b>SOLN STOP</b>	1	S-1039	2025-01	2027-01-31
<b>SUB TMB</b>	1	T-077	2022-09	2027-03-31
<b>DIL</b>	1	D-490	2024-12	2026-12-31

Component	Concentration	OD Value	OD Validation Criteria
<b>CAL A</b>	0 IU/mL	0.019	< 0.200
<b>CAL B</b>	50 IU/mL	0.907	> 0.300 & > Standard A
<b>CAL C</b>	100 IU/mL	1.470	> 0.500 & > Standard B
<b>CAL D</b>	200 IU/mL	2.330	> 1.000 & > Standard C

Final Report issued: 2025-04-25

Approved by: \_\_\_\_\_

Quality Control

Approved on: 2025-04-25



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany