



**GOLD
STANDARD
DIAGNOSTICS**



Certificate of Analysis (CoA)

Product: **Toxoplasma gondii IgG**
 Prod.No.: **RE57101**
 Lot: **TOXG-095**
 Expiry: **2027-08-31**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0460-095	2026-02	2027-08-31
CONJ	1	KG0460-095	2026-02	2027-08-31
CAL A	1	SAG0460-095	2026-02	2027-08-31
CAL B	1	SBG0460-095	2026-02	2027-08-31
CAL C	1	SCG0460-095	2026-02	2027-08-31
CAL D	1	SDG0460-095	2026-02	2027-08-31
WASH BUF 20x	1	W-1135	2026-02	2028-02-29
SOLN STOP	1	S-1067	2026-01	2028-01-31
SUB TMB	1	T-079	2025-12	2028-12-10
DIL	1	D-499	2026-01	2028-01-31

Component	Concentration	OD Value	OD Validation Criteria
CAL A	0 IU/mL	0.012	< 0.200
CAL B	50 IU/mL	0.571	> 0.300 & > Standard A
CAL C	100 IU/mL	0.987	> 0.500 & > Standard B
CAL D	200 IU/mL	1.679	> 1.000 & > Standard C

Final Report issued: 2026-03-04

Approved by: 
 Quality Control

Approved on: 2026-03-04



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany